

Recombinant Rat CD192/CCR2 Protein, N-GST & C-His

Catalog No: #AP98093

Package Size: #AP98093-1 50ug #AP98093-2 100ug #AP98093-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Rat CD192/CCR2 Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	C-C chemokine receptor type 2; C-C CKR-2; CC-CKR-2; CCR-2; CCR2; CD192; Ccr2; Cmkbr2
Accession No.	O55193
Target Species	Rattus norvegicus (Rat)
Calculated MW	34.90 kDa
Target Length	Met1-Pro60
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.