

Recombinant Mouse ORM1/a1AGP Protein, N-His

Catalog No: #AP98097

Package Size: #AP98097-1 50ug #AP98097-2 100ug #AP98097-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse ORM1/a1AGP Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	ORM1, Orosomucoid-1, AGP1, OMD 1, AGP 1, Alpha-1-acid glycoprotein 1,
Accession No.	Q60590
Target Species	Mus musculus (Mouse)
Calculated MW	24.23 kDa
Target Length	Gln19-Ala207
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.