

Recombinant BVDV NS3/p80 Protein, N-His

Catalog No: #AP98100

Package Size: #AP98100-1 50ug #AP98100-2 100ug #AP98100-3 1mg

Description

Product Name	Recombinant BVDV NS3/p80 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Genome polyprotein, Serine protease NS3, 3.4.21.113, 3.6.1.15, 3.6.4.13, Non-structural protein 3, p80
Accession No.	P19711
Target Species	Bovine viral diarrhea virus (isolate NADL) (BVDV) (Mucosal disease virus)
Calculated MW	18.81 kDa
Target Length	Glu1887-Thr2031
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.