

Recombinant BVDV NS2/p54 Protein, N-GST & C-His

Catalog No: #AP98101

Package Size: #AP98101-1 50ug #AP98101-2 100ug #AP98101-3 1mg

Description

Product Name	Recombinant BVDV NS2/p54 Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Genome polyprotein, Cysteine protease NS2, 3.4.22.-, Non-structural protein 2, p54
Accession No.	P19711
Target Species	Bovine viral diarrhea virus (isolate NADL) (BVDV) (Mucosal disease virus)
Calculated MW	37.55 kDa
Target Length	Met1537-Pro1617
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.