

Recombinant Human $\alpha 5\beta 1$ (ITGA5 & ITGB1) Protein,
C-Strep & C-His

Catalog No: #AP98111

Package Size: #AP98111-1 50ug #AP98111-2 100ug #AP98111-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human $\alpha 5\beta 1$ (ITGA5 & ITGB1) Protein, C-Strep & C-His
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	MSK12, VLA-4 subunit beta, CD29, Integrin beta-1, ITGB1, FNRB, MDF2, GPIIA, Glycoprotein IIa, Fibronectin receptor subunit beta, CD49 antigen-like family member E, Integrin alpha-F, ITGA5, VLA-5, Fibronectin receptor subunit alpha, CD49e, Integrin alpha-5, FNRA
Accession No.	P05556 & P08648
Target Species	Homo sapiens (Human)
Calculated MW	86.86 kDa & 109.56 kDa
Target Length	[Met1-Asp728] & [Met1-Tyr995]
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.