

Recombinant Human MMAB Protein, N-His

Catalog No: #AP98117

Package Size: #AP98117-1 50ug #AP98117-2 100ug #AP98117-3 1mg

Description

Product Name	Recombinant Human MMAB Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	MMAB, Methylmalonic aciduria type B protein, Cob(II)yrinic acid a,c-diamide adenosyltransferase, Cobinamide/cobalamin adenosyltransferase, Corrinoid adenosyltransferase, Cob(II)alamin adenosyltransferase
Accession No.	Q96EY8
Target Species	Homo sapiens (Human)
Calculated MW	23.16 kDa
Target Length	Met1-Ala186
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.