

Recombinant Dog CD137/4-1BB Protein, C-Fc

Catalog No: #AP98127

Package Size: #AP98127-1 50ug #AP98127-2 100ug #AP98127-3 1mg

Description

Product Name	Recombinant Dog CD137/4-1BB Protein, C-Fc
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Tumor necrosis factor receptor superfamily member 9, CD137, 4-1BB ligand receptor, CDw137, TNFRSF9, T-cell antigen ILA, T-cell antigen 4-1BB homolog, ILA,
Accession No.	A0A8C0SEC4
Target Species	Canis lupus familiaris (Dog) (Canis familiaris)
Calculated MW	63 kDa
Target Length	Ile24-Gln186
Tag Info	Fc-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.