

Recombinant Mycobacterium tuberculosis esat6 Protein, C-His

Catalog No: #AP98134

Package Size: #AP98134-1 50ug #AP98134-2 100ug #AP98134-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mycobacterium tuberculosis esat6 Protein, C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	6 kDa early secretory antigenic target, ESAT-6, esxA, esat6, Rv3875, MTV027.10
Accession No.	P9WNK7
Target Species	Mycobacterium tuberculosis (strain ATCC 25618 / H37Rv)
Calculated MW	10.73 kDa
Target Length	Met1-Ala95
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.