

Recombinant Human MRPL49 Protein, N-GST & C-His

Catalog No: #AP98144

Package Size: #AP98144-1 50ug #AP98144-2 100ug #AP98144-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|----------------|---|
| Product Name | Recombinant Human MRPL49 Protein, N-GST & C-His |
| Host Species | E. coli |
| Purification | >90% as determined by SDS-PAGE. |
| Applications | ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress |
| Other Names | L49mt, NOF1, MRPL49, Mitochondrial large ribosomal subunit protein mL49, 39S ribosomal protein L49, mitochondrial, C11orf4, NOF, Neighbor of FAU, Protein NOF1, MRP-L49 |
| Accession No. | Q13405 |
| Target Species | Homo sapiens (Human) |
| Calculated MW | 38.44 kDa |
| Target Length | Pro79-Phe166 |
| Formulation | Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol. |
| Storage | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.