

Recombinant Influenza A virus (H7N9) M2e Protein, N-GST & C-His

Catalog No: #AP98169

Package Size: #AP98169-1 50ug #AP98169-2 100ug #AP98169-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Influenza A virus (H7N9) M2e Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Matrix protein 2, Proton channel protein M2, M
Accession No.	AGI60299.2
Target Species	Influenza A virus (A/Hangzhou/1/2013(H7N9))
Calculated MW	30.90 kDa
Target Length	Met1-Pro25
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.