

Recombinant Influenza A virus (H5N1) NF (CsgA)-3M2e Protein, N-His

Catalog No: #AP98171

Package Size: #AP98171-1 50ug #AP98171-2 100ug #AP98171-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Influenza A virus (H5N1) NF (CsgA)-3M2e Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Matrix protein 2, Proton channel protein M2, M, Major curlin subunit, csgA, b1042, JW1025
Accession No.	P28307 & YP_308670.1
Target Species	Influenza A virus (strain A/Goose/Guangdong/1/1996 H5N1 genotype Gs/Gd)
Calculated MW	24.59 kDa
Target Length	[Gly21-Tyr151] & [Met1-Pro25]3
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.