

Recombinant IBV S1/Spike glycoprotein 1 (NTD) Protein, C-Fc

Catalog No: #AP98182

Package Size: #AP98182-1 50ug #AP98182-2 100ug #AP98182-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant IBV S1/Spike glycoprotein 1 (NTD) Protein, C-Fc
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Spike glycoprotein, S glycoprotein, E2, Peplomer protein, Spike protein S1
Accession No.	P12651
Target Species	Avian infectious bronchitis virus (strain M41) (IBV)
Calculated MW	51.76 kDa
Target Length	Tyr21-Asn237
Tag Info	Fc-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.