

## Recombinant Tupaia chinensis CD69 Protein, C-His

Catalog No: #AP98209

Package Size: #AP98209-1 50ug #AP98209-2 100ug #AP98209-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Tupaia chinensis CD69 Protein, C-His
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	AIM, CD69, C-type lectin domain family 2 member C, Early T-cell activation antigen p60, BL-AC/P26, GP32/28, MLR-3, Early activation antigen CD69, Leukocyte surface antigen Leu-23, CLEC2C, EA1, Activation inducer molecule,
Accession No.	XP_006141645.1
Target Species	Tupaia chinensis (Chinese tree shrew)
Calculated MW	18.52 kDa
Target Length	Ser62-Gln199
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.