

Recombinant Tupaia chinensis CD69 Protein, C-Fc

Catalog No: #AP98210

Package Size: #AP98210-1 50ug #AP98210-2 100ug #AP98210-3 1mg

Description

Product Name	Recombinant Tupaia chinensis CD69 Protein, C-Fc
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	AIM, CD69, C-type lectin domain family 2 member C, Early T-cell activation antigen p60, BL-AC/P26, GP32/28, MLR-3, Early activation antigen CD69, Leukocyte surface antigen Leu-23, CLEC2C, EA1, Activation inducer molecule,
Accession No.	XP_006141645.1
Target Species	Tupaia chinensis (Chinese tree shrew)
Calculated MW	43.95 kDa
Target Length	Ser62-Gln199
Tag Info	Fc-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.