

Recombinant Tupaia chinensis CD152/CTLA4 Protein, C-Strep

Catalog No: #AP98215

Package Size: #AP98215-1 50ug #AP98215-2 100ug #AP98215-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Tupaia chinensis CD152/CTLA4 Protein, C-Strep
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	CD152, Cytotoxic T-lymphocyte protein 4, CTLA-4, CTLA4, Cytotoxic T-lymphocyte-associated antigen 4,
Accession No.	XP_006153140.1
Target Species	Tupaia chinensis (Chinese tree shrew)
Calculated MW	16.48 kDa
Target Length	Met38-Asp161
Tag Info	Strep-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.