

Recombinant Dog OSMR/IL-31RB Protein, C-hFc

Catalog No: #AP98228

Package Size: #AP98228-1 50ug #AP98228-2 100ug #AP98228-3 1mg

Description

Product Name	Recombinant Dog OSMR/IL-31RB Protein, C-hFc
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	OSMRB, IL-31R subunit beta, Oncostatin-M-specific receptor subunit beta, Interleukin-31 receptor subunit beta, IL-31RB, IL-31R-beta, OSMR, IL-31 receptor subunit beta,
Accession No.	XP_005619431.1
Target Species	Canis lupus familiaris (Dog) (Canis familiaris)
Calculated MW	109.08 kDa
Target Length	Glu24-Glu732
Tag Info	Fc-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.