

Recombinant Human FSH Dimer Protein, C-His

Catalog No: #AP98232

Package Size: #AP98232-1 50ug #AP98232-2 100ug #AP98232-3 1mg

Description

Product Name	Recombinant Human FSH Dimer Protein, C-His
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Follicle-stimulating hormone alpha chain, Lutropin alpha chain, Thyrotropin alpha chain, Folitropin alpha chain, Luteinizing hormone alpha chain, Thyroid-stimulating hormone alpha chain, FSH-alpha, LSH-alpha, Glycoprotein hormones alpha chain, CGA, TSH-alpha, Choriogonadotropin alpha chain, Anterior pituitary glycoprotein hormones common subunit alpha, CG-alpha, Chorionic gonadotrophin subunit alpha
Accession No.	P01215 & P01225
Target Species	Homo sapiens (Human)
Calculated MW	20.64 kDa
Target Length	[Met1-Ser116] & [Met1-Glu129]
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.