

Recombinant Honeybee Api m 11/MRJP1 Protein, N-His-SUMO

Catalog No: #AP98253

Package Size: #AP98253-1 50ug #AP98253-2 100ug #AP98253-3 1mg

Description

Product Name	Recombinant Honeybee Api m 11/MRJP1 Protein, N-His-SUMO
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Major royal jelly protein 1, MRJP-1, MRJP1, 56-kDa protein 4, p56kP-4, Apalbumin 1, Apisin subunit MRJP1, Bee-milk protein, Royalactin, Jellein-1, Jelleine-I, Jellein-2, Jelleine-II, Jellein-4, Jelleine-IV, MRJP1
Accession No.	O18330
Target Species	Apis mellifera (Honeybee)
Calculated MW	59.39 kDa
Target Length	Asn20-Leu432
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.