

Recombinant Influenza A virus (H1N1) HA/Hemagglutinin Protein, C-His

Catalog No: #AP98311

Package Size: #AP98311-1 50ug #AP98311-2 100ug #AP98311-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Influenza A virus (H1N1) HA/Hemagglutinin Protein, C-His
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Hemagglutinin, Hemagglutinin HA1 chain, Hemagglutinin HA2 chain, HA
Accession No.	P03452
Target Species	Influenza A virus (strain A/Puerto Rico/8/1934 (H1N1))
Calculated MW	64.36 kDa
Target Length	Asp18-Gly518
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.