

# Recombinant Influenza A virus (H1N1) HA/Hemagglutinin Protein, C-8His

Catalog No: #AP98349

Package Size: #AP98349-1 50ug #AP98349-2 100ug #AP98349-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                |                                                                                                                                                                       |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name   | Recombinant Influenza A virus (H1N1) HA/Hemagglutinin Protein, C-8His                                                                                                 |
| Host Species   | Mammalian cells                                                                                                                                                       |
| Purification   | >90% as determined by SDS-PAGE.                                                                                                                                       |
| Applications   | ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress                                                                                                       |
| Other Names    | Hemagglutinin, Hemagglutinin HA1 chain, Hemagglutinin HA2 chain, HA                                                                                                   |
| Accession No.  | WHO42457.1                                                                                                                                                            |
| Target Species | Influenza A virus (A/Utah/11/2023 (H1N1))                                                                                                                             |
| Calculated MW  | 65.10 kDa                                                                                                                                                             |
| Target Length  | Ala17-Ile530                                                                                                                                                          |
| Formulation    | Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.                                                                                       |
| Storage        | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.