

Recombinant Influenza A virus (H5N6) HA/Hemagglutinin Protein, C-6His

Catalog No: #AP98372

Package Size: #AP98372-1 50ug #AP98372-2 100ug #AP98372-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Influenza A virus (H5N6) HA/Hemagglutinin Protein, C-6His
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Hemagglutinin, Hemagglutinin HA1 chain, Hemagglutinin HA2 chain, HA
Accession No.	A0A1Z4EAT9
Target Species	Influenza A virus (A/tundra swan/Tottori/3111S002/2016 (H5N6))
Calculated MW	64.23 kDa
Target Length	Met1-Ile530
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.