

# Recombinant Influenza A virus (H5N6) NA/Neuraminidase Protein, C-His

Catalog No: #AP98373

Package Size: #AP98373-1 50ug #AP98373-2 100ug #AP98373-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                |                                                                                                                                                                       |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name   | Recombinant Influenza A virus (H5N6) NA/Neuraminidase Protein, C-His                                                                                                  |
| Host Species   | Mammalian cells                                                                                                                                                       |
| Purification   | >90% as determined by SDS-PAGE.                                                                                                                                       |
| Applications   | ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress                                                                                                       |
| Other Names    | neuraminidase, NA                                                                                                                                                     |
| Accession No.  | AIK26468.1                                                                                                                                                            |
| Target Species | Influenza A virus (A/Anas_crecca/Hubei/Chenhu1623-5/2014 (H5N6))                                                                                                      |
| Calculated MW  | 51.67 kDa                                                                                                                                                             |
| Target Length  | Gly30-Lys470                                                                                                                                                          |
| Tag Info       | His-tag                                                                                                                                                               |
| Formulation    | Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.                                                                                       |
| Storage        | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.