

Recombinant Human NPY1R Protein, N-GST & C-His

Catalog No: #AP98395

Package Size: #AP98395-1 50ug #AP98395-2 100ug #AP98395-3 1mg

Description

Product Name	Recombinant Human NPY1R Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	NPYY1, NPYR, NPY1R, NPY1-R, Neuropeptide Y receptor type 1
Accession No.	P25929
Target Species	Homo sapiens (Human)
Calculated MW	37.49 kDa
Target Length	[Asn2-Ala44] & [Tyr176-Tyr211]
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.