

D-Lactyl Lysine Rabbit mAb

Catalog No: #HL002

Package Size: #HL002-1 100ul



Description

Product Name	D-Lactyl Lysine Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purification	>95%. Protein A/G purified from cell culture supernatant.
Applications	WB;IF;ChIP
Species Reactivity	Species independent
Conjugates	Unconjugated
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20Λ C

Application Details

WB 1:500-1:2000; IHC 1:100-1:300; ChIP 1:20-200;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.