

Histone H2A.Z (L-Lactyl-Lys11) Rabbit mAb

Catalog No: #HL003

Package Size: #HL003-1 100ul

Description

Product Name	Histone H2A.Z (L-Lactyl-Lys11) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purification	>95%. Protein A/G purified from cell culture supernatant.
Applications	WB;IHC;IF;IP
Species Reactivity	Human;Mouse;Rat
Conjugates	Unconjugated
Target Name	H2AZ1
Other Names	H2A.Z; H2AFZ; MGC117173
Accession No.	P0C0S5
Calculated MW	14 kDa
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20Λ C

Application Details

WB 1:500-1:2000; IHC 1:100-1:300; IF 1:100-1:300; IP 1:20-500;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.