

Goat Anti-Chicken IgY (H&L) HRP

Catalog No: #L35001

Package Size: #L35001 100 ul



Description

Product Name	Goat Anti-Chicken IgY (H&L) HRP
Host Species	Goat
Clonality	Polyclonal
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, WB, IH, IC
Species Reactivity	C
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Chicken IgY. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Chicken IgY
Conjugates	HRP
Concentration	1 mg/ml
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide

Application Details

E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)

Published Papers

el at., SynCAM1 deficiency in the hippocampal parvalbumin interneurons contributes to sevoflurane-induced cognitive impairment in neonatal ratsIn CNS Neurosci Ther On2023 Dec 17byMing-Ming Zhao , Han-Wen Gu et al..PMID:38105652, , (2023)
[PMID:38105652](#)

el at., SynCAM1 deficiency in the hippocampal parvalbumin interneurons contributes to sevoflurane-induced cognitive impairment in neonatal rats. In CNS Neurosci Ther on 2024 Jan by Ming-Ming Zhao, Han-Wen Gu,et al..PMID:38105652, , (2024)
[PMID:38105652](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.